

TABLETS & CAPSULES

www.tabletsandcapsules.com

May 2018

Volume 16 Number 4

Cover photo shows a multi-product track feeder on a Hoonga HMV3 blister packaging machine feeding several different capsules into a single blister pack. Feeding



options include a brush box or shutter gate. Company is the sole North American distributor and service provider of the Hoonga line of flexible, economic blister machines that are manufactured in Hoonga's 40,000-square-foot factory in Ojeong, Korea. Units handle tablets, capsules, and softgels and outputs range from 150 to 900 blisters per minute. Photo courtesy of Key International, Cranbury, NJ. For more information, call 609 619 3685 or visit www.keyinternational.com.

packaging

From serialization to aggregation: A holistic approach

8

Jörg Willburger

Serialization ready?

12

Nate Todd

Issue focus: Packaging

18

also in this issue

Using vacuum transfer equipment to improve pharmaceutical processing

21

Doan Pendleton

Powder transfer & containment: A survey of what's available

25

How to get better grinding performance with pneumatic conveying

28

Gary Liu

Particle size reduction separation & analysis: A survey of what's available

33

CRS annual meeting preview

36

Eye on Excipients

38

Hanna Bogner and Gernot Warnke

Back page

Cover 3

Wesley Mancoff

departments

Research

4

People news

43

Markets

4

Supplier news

44

Industry news

4

Product update

44

Calendar

42

Advertiser index

44

Call for papers

43